

NEBNext[®] Multiplex Oligos for Illumina[®] (Index Primers Set 2)

NEB #E7500S/L

24/96 reactions

Version 5.0_1/20

Table of Contents

| | |
|---|---|
| NEBNext Adaptors and Primers for Illumina | 3 |
| Kit Components | 4 |
| Revision History | 5 |

The NEBNext Multiplex Oligos for Illumina (Index Primers Set 2) Includes:

The volumes provided are sufficient for preparation of up to 24 reactions (NEB #E7500S) and 96 reactions (NEB #E7500L). All reagents should be stored at –20°C.

NEBNext Adaptor for Illumina

USER[®] Enzyme

NEBNext Universal PCR Primer for Illumina

NEBNext Index 13 Primer for Illumina

NEBNext Index 14 Primer for Illumina

NEBNext Index 15 Primer for Illumina

NEBNext Index 16 Primer for Illumina

NEBNext Index 18 Primer for Illumina

NEBNext Index 19 Primer for Illumina

NEBNext Index 20 Primer for Illumina

NEBNext Index 21 Primer for Illumina

NEBNext Index 22 Primer for Illumina

NEBNext Index 23 Primer for Illumina

NEBNext Index 25 Primer for Illumina

NEBNext Index 27 Primer for Illumina

For the list of additional materials required, please check the manual for your NEBNext Library Prep Kit.

Applications


The NEBNext Multiplex Oligos for Illumina (Index Primers Set 2) contains adaptors and primers that are ideally suited for multiplex sample preparation for next generation sequencing on the Illumina platform (Illumina, Inc.). Each of these components must pass rigorous quality control standards and is lot controlled, both individually and as a set of reagents.

Lot Control: The lots provided in the NEBNext Multiplex Oligos for Illumina (Index Primers Set 2) are managed separately and are qualified by additional functional validation. Individual reagents undergo standard enzyme activity and quality control assays, and also meet stringent criteria in the additional quality controls listed on each individual component page.

Functionally Validated: Each of the components is functionally validated through construction and sequencing of a genomic DNA libraries on the Illumina platform.

For larger volume requirements, customized and bulk packaging is available by purchasing through the OEM/Bulks department at NEB. Please contact OEM@neb.com for further information.

Figure 1. Workflow



Please Refer to the Kit Specific Protocol for using the NEBNext Multiplex Oligos for Illumina

The following kits are designed for use with the NEBNext Multiplex Oligos for Illumina:

- #E7760, NEBNext Ultra™ II Directional RNA Library Prep Kit for Illumina
- #E7765, NEBNext Ultra II Directional RNA Library Prep with Sample Purification Beads
- #E7770, NEBNext Ultra II RNA Library Prep Kit for Illumina
- #E7775, NEBNext Ultra II RNA Library Prep with Sample Purification Beads
- #E7805, NEBNext Ultra II FS DNA Library Prep Kit for Illumina
- #E6177, NEBNext Ultra II FS DNA Library Prep with Sample Purification Beads
- #E7645, NEBNext Ultra II DNA Library Prep Kit for Illumina
- #E7103, NEBNext Ultra II DNA Library Prep with Sample Purification Beads
- #E7595, NEBNext Ultra II Ligation Module
- #E7420, NEBNext Ultra Directional RNA Library Prep Kit for Illumina
- #E7530, NEBNext Ultra RNA Library Prep Kit for Illumina
- #E7370, NEBNext Ultra DNA Library Prep Kit for Illumina
- #E7445, NEBNext Ultra Ligation Module
- #E6040, NEBNext DNA Library Prep Master Mix Set for Illumina
- #E6240, NEBNext ChIP-Seq Library Prep Master Mix Set for Illumina
- #E6056, NEBNext Quick Ligation Module

NEBNext Adaptors and Primers for Illumina

| INDEX PRIMER | INDEX PRIMER SEQUENCE | EXPECTED INDEX SEQUENCE READ |
|---|---|------------------------------|
| NEBNext Index 13 Primer for Illumina | 5'-CAAGCAGAAGACGGCATAACGAGATTG <u>TTGACT</u> GTGACTGGAGTTCAGACGTGTGCTCTCCGATC-s-T-3' | AGTCAA |
| NEBNext Index 14 Primer for Illumina | 5'-CAAGCAGAAGACGGCATAACGAGATAC <u>GGAAC</u> TGTGACTGGAGTTCAGACGTGTGCTCTCCGATC-s-T-3' | AGTTCC |
| NEBNext Index 15 Primer for Illumina | 5'-CAAGCAGAAGACGGCATAACGAGATT <u>CTGACAT</u> GTGACTGGAGTTCAGACGTGTGCTCTCCGATC-s-T-3' | ATGTCA |
| NEBNext Index 16 Primer for Illumina | 5'-CAAGCAGAAGACGGCATAACGAGATGC <u>GGACGG</u> TGTGACTGGAGTTCAGACGTGTGCTCTCCGATC-s-T-3' | CCGTCC |
| NEBNext Index 18 Primer for Illumina | 5'-CAAGCAGAAGACGGCATAACGAGATGT <u>GCGGAC</u> TGTGACTGGAGTTCAGACGTGTGCTCTCCGATC-s-T-3' | GTCCGC |
| NEBNext Index 19 Primer for Illumina | 5'-CAAGCAGAAGACGGCATAACGAGATCG <u>TTTCAC</u> TGTGACTGGAGTTCAGACGTGTGCTCTCCGATC-s-T-3' | GTGAAA |
| NEBNext Index 20 Primer for Illumina | 5'-CAAGCAGAAGACGGCATAACGAGATAA <u>GCCAC</u> TGTGACTGGAGTTCAGACGTGTGCTCTCCGATC-s-T-3' | GTGGCC |
| NEBNext Index 21 Primer for Illumina | 5'-CAAGCAGAAGACGGCATAACGAGATT <u>CGAAAC</u> TGTGACTGGAGTTCAGACGTGTGCTCTCCGATC-s-T-3' | GTTTCG |
| NEBNext Index 22 Primer for Illumina | 5'-CAAGCAGAAGACGGCATAACGAGATT <u>AGTACG</u> TGTGACTGGAGTTCAGACGTGTGCTCTCCGATC-s-T-3' | CGTACG |
| NEBNext Index 23 Primer for Illumina | 5'-CAAGCAGAAGACGGCATAACGAGATAT <u>CCACTC</u> TGTGACTGGAGTTCAGACGTGTGCTCTCCGATC-s-T-3' | GAGTGG |
| NEBNext Index 25 Primer for Illumina | 5'-CAAGCAGAAGACGGCATAACGAGATAT <u>ATCAGT</u> TGTGACTGGAGTTCAGACGTGTGCTCTCCGATC-s-T-3' | ACTGAT |
| NEBNext Index 27 Primer for Illumina | 5'-CAAGCAGAAGACGGCATAACGAGATAA <u>AGGAAT</u> TGTGACTGGAGTTCAGACGTGTGCTCTCCGATC-s-T-3' | ATTCTT |
| NEBNext Adaptor for Illumina | 5'-/5Phos/GAT CGG AAG AGC ACA CGT CTG AAC TCC AGT CUA CAC TCT TTC CCT ACA CGA CGC TCT TCC GAT C-s-T-3' | N/A |
| NEBNext Universal PCR Primer for Illumina | 5'-AAT GAT ACG GCG ACC ACC GAG ATC TAC ACT CTT TCC CTA CAC GAC GCT CTT CCG ATC-s-T-3' | N/A |

Where -s- indicates phosphorothioate bond.

Note: If fewer than 12 indexes are used in a lane for sequencing, it is recommended to use the following indexes:

Pool of 2 samples: Index #18 and 25

Pool of 3 samples: Index #13, 18 and 23

Pool of 4 samples: Index #16, 18, 20 and 23

Kit Components

NEB #E7500S Table of Components

| NEB # | PRODUCT | VOLUME |
|--------|---|----------|
| E7337A | NEBNext Adaptor for Illumina | 0.24 ml |
| E7338A | USER Enzyme | 0.072 ml |
| E6861A | NEBNext Universal PCR Primer for Illumina | 0.120 ml |
| E6862A | NEBNext Index 13 Primer for Illumina | 0.010 ml |
| E6863A | NEBNext Index 14 Primer for Illumina | 0.010 ml |
| E6864A | NEBNext Index 15 Primer for Illumina | 0.010 ml |
| E7511A | NEBNext Index 16 Primer for Illumina | 0.010 ml |
| E6866A | NEBNext Index 18 Primer for Illumina | 0.010 ml |
| E6867A | NEBNext Index 19 Primer for Illumina | 0.010 ml |
| E6868A | NEBNext Index 20 Primer for Illumina | 0.010 ml |
| E6869A | NEBNext Index 21 Primer for Illumina | 0.010 ml |
| E6870A | NEBNext Index 22 Primer for Illumina | 0.010 ml |
| E6871A | NEBNext Index 23 Primer for Illumina | 0.010 ml |
| E6872A | NEBNext Index 25 Primer for Illumina | 0.010 ml |
| E6873A | NEBNext Index 27 Primer for Illumina | 0.010 ml |

NEB #E7500L Table of Components

| NEB # | PRODUCT | VOLUME |
|---------|---|----------|
| E7337AA | NEBNext Adaptor for Illumina | 0.96 ml |
| E7338AA | USER Enzyme | 0.288 ml |
| E6861AA | NEBNext Universal PCR Primer for Illumina | 0.48 ml |
| E6862AA | NEBNext Index 13 Primer for Illumina | 0.04 ml |
| E6863AA | NEBNext Index 14 Primer for Illumina | 0.04 ml |
| E6864AA | NEBNext Index 15 Primer for Illumina | 0.04 ml |
| E7511AA | NEBNext Index 16 Primer for Illumina | 0.04 ml |
| E6866AA | NEBNext Index 18 Primer for Illumina | 0.04 ml |
| E6867AA | NEBNext Index 19 Primer for Illumina | 0.04 ml |
| E6868AA | NEBNext Index 20 Primer for Illumina | 0.04 ml |
| E6869AA | NEBNext Index 21 Primer for Illumina | 0.04 ml |
| E6870AA | NEBNext Index 22 Primer for Illumina | 0.04 ml |
| E6871AA | NEBNext Index 23 Primer for Illumina | 0.04 ml |
| E6872AA | NEBNext Index 25 Primer for Illumina | 0.04 ml |
| E6873AA | NEBNext Index 27 Primer for Illumina | 0.04 ml |

Revision History

| REVISION # | DESCRIPTION | DATE |
|------------|--|-------|
| 1.0 | N/A | 1/13 |
| 2.0 | Change Index Primer and Universal Primer concentration from 25 μ M to 10 μ M. All part #'s have changed. Concentrations of Universal PCR Primer and Index Primers have changed from 25 μ M to 10 μ M. | 1/15 |
| 3.0 | Switched out Index 16 Primer with a new sequence. | 10/15 |
| 3.1 | Update "Required Materials Not Included for DNA or ChIP Libraries" and "Required Materials Not Included for RNA Libraries". Update the list of kits using the NEBNext Multiplex Oligos for Illumina. | 4/16 |
| 4.0 | Delete "Required Materials Not Included for DNA or ChIP Libraries" and "Required Materials Not Included for RNA Libraries" and replace with guidance note. Update the list of kits using the NEBNext Multiplex Oligos for Illumina. Note, Pool of 4 Samples: Index# 13, 14, 16 and 18, replace with 16, 18, 20 and 23. DNA sequence for NEBNext Adaptor for Illumina has been corrected. Change part number for NEBNext Index 16 Primer for Illumina and add 2 additional bases to the sequence. Create "Kit Component – Table of Components" for small and large size kits. Change Index Primer sequences table to include adaptor, and NEBNext Universal PCR Primer and rename it as NEBNext Adaptors and Primers for Illumina. Delete individual component information pages. Updated the functionally validated statement on page 2. | 1/18 |
| 5.0 | Update to new manual format | 1/20 |

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

This product is covered by one or more patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information about commercial rights, please email us at ebd@neb.com. While NEB develops and validates its products for various applications, the use of this product may require the buyer to obtain additional third party intellectual property rights for certain applications.

AMPURE® and SPRISELECT® are registered trademarks of Beckman Coulter, Inc.

BIOANALYZER® is a registered trademark of Agilent Technologies, Inc.

EPPENDORF® and LOBIND® are registered trademarks of Eppendorf AG.

ILLUMINA® is a registered trademark of Illumina, Inc.

E-GEL® is a registered trademark of Life Technologies, Inc.

© Copyright 2020, New England Biolabs, Inc.; all rights reserved



www.neb.com



be INSPIRED
drive DISCOVERY
stay GENUINE

New England Biolabs, Inc., 240 County Road, Ipswich, MA 01938-2723 Telephone: (978) 927-5054 Toll Free: (USA Orders) 1-800-632-5227 (USA Tech) 1-800-632-7799 Fax: (978) 921-1350 e-mail: info@neb.com