240 County Road Ipswich, MA 01938-2723 Tel 978-927-5054 Fax 978-921-1350

Date

30 Mar 2016

www.neb.com info@neb.com

New England Biolabs Product Specification

Product Name: Lemo21(DE3) Competent E. coli

Catalog #: C2528H Shelf Life: 12 months Storage Temp: -80°C

Specification Version: PS-C2528H v1.0 Effective Date: 30 Mar 2016

Assay Name/Specification (minimum release criteria)

Antibiotic Resistance (Chloramphenicol) - 15 µl of untransformed Lemo21(DE3) Competent E. coli streaked onto a Rich Broth plate containing Chloramphenicol will form colonies after incubation for 16 hours at 37°C.

Antibiotic Sensitivity (Ampicillin) - 15 µl of untransformed Lemo21(DE3) Competent E. coli streaked onto a Rich Broth plate containing Ampicillin will not form colonies after incubation for 16 hours at 37°C.

Antibiotic Sensitivity (Kanamycin) - 15 µl of untransformed Lemo21(DE3) Competent E. coli streaked onto a Rich Broth plate containing Kanamycin will not form colonies after incubation for 16 hours at 37°C.

Antibiotic Sensitivity (Spectinomycin) - 15 µl of untransformed Lemo21(DE3) Competent E. coli streaked onto a Rich Broth plate containing Spectinomycin will not form colonies after incubation for 16 hours at 37°C.

Antibiotic Sensitivity (Streptomycin) - 15 μl of untransformed Lemo21(DE3) Competent *E. coli* streaked onto a Rich Broth plate containing Streptomycin will not form colonies after incubation for 16 hours at 37°C.

Antibiotic Sensitivity (Tetracycline) - 15 µl of untransformed Lemo21(DE3) Competent E. coli streaked onto a Rich Broth plate containing Tetracycline will not form colonies after incubation for 16 hours at 37°C.

Phage Resistance (Ф 80) - 15 µl of untransformed Lemo21(DE3) Competent E. coli streaked onto a Rich Broth plate does not support plaque formation by phage Φ 80 after incubation for 16 hours at 37°C.

Transformation Efficiency - 50 µl of Lemo21(DE3) Competent E. coli cells were transformed with 100 pg of pUC19 DNA using the transformation protocol provided. Incubation overnight on LB-Ampicillin plates at 37°C resulted in >1 x 10e7 cfu/µg of DNA.

Derek Robinson

Director of Quality Control





