

240 County Road Ipswich, MA 01938-2723 Tel 978-927-5054 Fax 978-921-1350 www.neb.com info@neb.com

New England Biolabs Product Specification

Product Name:	NiCo21(DE3) Competent E. coli
Catalog #:	С2529Н
Shelf Life:	12 months
Storage Temp:	-80°C
Specification Version:	PS-C2529H v1.0
Effective Date:	30 Mar 2016

Assay Name/Specification (minimum release criteria)

Antibiotic Sensitivity (Ampicillin) - 15 µl of untransformed NiCo21(DE3) Competent *E. coli* streaked onto a Rich Broth plate containing Ampicillin will not form colonies after incubation for 16 hours at 37°C.

Antibiotic Sensitivity (Chloramphenicol) - 15 µl of untransformed NiCo21(DE3) Competent *E. coli* streaked onto a Rich Broth plate containing Chloramphenicol will not form colonies after incubation for 16 hours at 37°C.

Antibiotic Sensitivity (Kanamycin) - 15 μ l of untransformed NiCo21(DE3) Competent *E. coli* streaked onto a Rich Broth plate containing Kanamycin will not form colonies after incubation for 16 hours at 37°C.

Antibiotic Sensitivity (Spectinomycin) - 15 μ l of untransformed NiCo21(DE3) Competent *E. coli* streaked onto a Rich Broth plate containing Spectinomycin will not form colonies after incubation for 16 hours at 37°C.

Antibiotic Sensitivity (Streptomycin) - 15 µl of untransformed NiCo21(DE3) Competent *E. coli* streaked onto a Rich Broth plate containing Streptomycin will not form colonies after incubation for 16 hours at 37°C.

Antibiotic Sensitivity (Tetracycline) - 15 μ l of untransformed NiCo21(DE3) Competent *E. coli* streaked onto a Rich Broth plate containing Tetracycline will not form colonies after incubation for 16 hours at 37°C.

Phage Resistance (Φ 80) - 15 µl of untransformed NiCo21(DE3) Competent *E. coli* streaked onto a Rich Broth plate does not support plaque formation by phage Φ 80 after incubation for 16 hours at 37°C.

Transformation Efficiency - 50 μ l of NiCo21(DE3) Competent *E. coli* cells were transformed with 100 pg of pUC19 DNA using the transformation protocol provided. Incubation overnight on LB-Ampicillin plates at 37°C resulted in >1 x 10e7 cfu/ μ g of DNA.

Date 30 Mar 2016

Derek Robinson Director of Quality Control



PS-C2529H v1.0 Page 1 of 1