

E. coli ER2420/ pACYC177



1-800-632-7799
info@neb.com
www.neb.com



E4151S 008160621061

E4151S

Lot: 0081606

Exp: 6/21

Store at -20°C

Description: A suspension of *E. coli* strain ER2420 containing the plasmid pACYC177 which has been grown in rich medium and brought to 50% glycerol. ER2420 background: *EcoK* r⁻ m⁻ / McrBC⁻ Mrr Dcm⁺.

Genotype: F⁻ *ara-14 leu fhuA2 Δ(gpt-proA)62 lacY1 glnV44 galK2 rpsL20(Str^R) xyl-5 mtl-1 Δ(mcrC-mrr)*_{HB101}

Antibiotic Resistance: The plasmid pACYC177 carries a gene encoding resistance to ampicillin (beta-lactamase) and a gene encoding resistance to kanamycin (kanamycin phosphotransferase). To maintain the plasmid cells should be grown with 50–100 µg/ml of amp and 35–50 µg/ml of kanamycin. The plasmid, which is present at a low copy number, should be amplified using a standard chloramphenicol procedure.

Notes: Storage at -70°C is recommended for periods longer than 30 days. Avoid repeated freeze/thaw cycles.

References:

1. Chang, A. and Cohen, S. (1978) *J. Bacteriol.* 134, 1141–1156.
2. Rose, R.E. (1988) *Nucleic Acids. Res.*16, 355.



NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

CERTIFICATE OF ANALYSIS

E. coli ER2420/ pACYC177



1-800-632-7799
info@neb.com
www.neb.com



E4151S 008160621061

E4151S

Lot: 0081606

Exp: 6/21

Store at -20°C

Description: A suspension of *E. coli* strain ER2420 containing the plasmid pACYC177 which has been grown in rich medium and brought to 50% glycerol. ER2420 background: *EcoK* r⁻ m⁻ / McrBC⁻ Mrr Dcm⁺.

Genotype: F⁻ *ara-14 leu fhuA2 Δ(gpt-proA)62 lacY1 glnV44 galK2 rpsL20(Str^R) xyl-5 mtl-1 Δ(mcrC-mrr)*_{HB101}

Antibiotic Resistance: The plasmid pACYC177 carries a gene encoding resistance to ampicillin (beta-lactamase) and a gene encoding resistance to kanamycin (kanamycin phosphotransferase). To maintain the plasmid cells should be grown with 50–100 µg/ml of amp and 35–50 µg/ml of kanamycin. The plasmid, which is present at a low copy number, should be amplified using a standard chloramphenicol procedure.

Notes: Storage at -70°C is recommended for periods longer than 30 days. Avoid repeated freeze/thaw cycles.

References:

1. Chang, A. and Cohen, S. (1978) *J. Bacteriol.* 134, 1141–1156.
2. Rose, R.E. (1988) *Nucleic Acids. Res.*16, 355.



NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

CERTIFICATE OF ANALYSIS