

50 ml Magnetic Separation Rack



1-800-632-7799
info@neb.com
www.neb.com



S1507S

holds 4-tubes

Lot: 0021209

Description: The 50 ml Magnetic Separation Rack is designed to be used for small-scale separations using magnetic particles.

Capacity: four 50 ml tubes.

Magnets: Neodymium rare earth permanent magnets.

CAUTION: Rare earth magnets can be extremely powerful. Care should be taken when handling these magnets. Keep magnetized parts away from mechanical/electrical instruments. The strong magnetic fields of neodymium magnets can also damage magnetic media such as floppy disks, credit cards, magnetic I.D. cards, cassette tapes, video tapes or other such devices. They can also damage televisions, VCRs, computer monitors and other CRT displays. Never place neodymium magnets near electronic appliances.

Never allow neodymium magnets near a person with a pacemaker or similar medical aid. The strong magnetic fields of the magnet can affect the operation of such devices.

Neodymium magnets will lose their magnetic properties if heated above 80°C (175°F). Neodymium magnets should never be burned, as burning them will create toxic fumes.



NEB's Magnetic Separation Racks are now available in three sizes: 6 or 12 (1.5 ml) tubes and 4 (50 ml) tubes. To order please see companion products listed below.

Companion Products:

Streptavidin Magnetic Beads #S1420S	20 mg
Protein A Magnetic Beads #S1425S	1 ml
Protein G Magnetic Beads #S1430S	1 ml
Goat Anti-Mouse IgG Magnetic Beads #S1431S	20 mg
Goat Anti-Rabbit IgG Magnetic Beads #S1432S	20 mg
Amylose Magnetic Beads #E8035S	25 mg
Chitin Magnetic Beads #E8036S	2 ml
6-Tube Magnetic Separation Rack #S1506S	holds 6-tubes
12-Tube Magnetic Separation Rack #S1509S	holds 12-tubes